

# ZEPHIR™ CPAP FACE MASK



Mask with safety valve



Mask without safety valve

- **Latex - Free**
- **Available in 3 sizes**
- **PVC Material**

## Features and Benefits

- **Face mask devised for CPAP, with valved inspiratory and expiratory ports.**
- **Thanks to the inflatable cushion, the mask adjusts perfectly to the patient's face preventing air leakage.**
- **High transparency allows good visibility.**
- **Comes with or without the safety valve.**
- **Expiratory port fits KOO adjustable PEEP valve, part #KM-809**
- **Single patient use.**




KOO MEDICAL EQUIPMENT

# ZEPHIR™ CPAP FACE MASK

## Technical Characteristics

Product Number	Product Description	Mask Size	Quantity Per Box
KM-301	Dual port CPAP Face Mask	Large	22
KM-301D	Dual port CPAP Face Mask	Large Deep	22
KM-303	Dual port CPAP Face Mask	Small	24
KM-301A	Dual port safety CPAP Face Mask	Large	22
KM-301DA	Dual port safety CPAP Face Mask	Large Deep	22
KM-303A	Dual port safety CPAP Face Mask	Small	24



**NOTE :**  
- Single patient use

Distributed by:



*Latex free*



**KOO MEDICAL EQUIPMENT**

**KOO EUROPE SRL**  
Via della Praia 12/D  
10090 - Buttigliera Alta (TORINO)  
ITALY  
Tel : +39 011 934 8192  
Fax : +39 011 934 8627  
Info@kooeurope.com

**KOO AMERICAS, INC**  
301 East Wallace Kneeland Blvd  
Building 224 - 112  
Shelton, WA 98584  
Tel : +360 426 3669  
Fax : +360 358 1338  
andrea@koo-shanghai.com

**KOO ASIA LTD.**  
Unit 1, 16/F, Cable TV Tower,  
9 Hoi Shing Road, Tsuen Wan,  
N.T. HONG KONG  
Tel : +852 2110 4780  
Fax : +852 2110 4783  
sales@koo-asia.com